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1. APPELBAUM P.C., LEATHERS D.J.  
Evaluation of the Rapid NFT System for Identification of Gram-Negative, Nonfermenting Rods. (1984) J. Clin. Microbiol. 20, 730-734.
2. BERNARDS A.T., VAN DER TOORN J., VAN BOVEN C.P.A., DIJKSHOORN L.  
Evaluation of the Ability of a Commercial System to Identify *Acinetobacter* Genomic Species. (1996) Eur. J. Clin. Microbiol. Infect. Dis. 15, 4, 303-308.
3. BILKEY M.K., BREMNER D.A., CAMERON G.L., GARNER J.G.  
Comparison of Five Commercial Methods for the Identification of Non-fermentative and Oxidase Positive Fermentative Gram-Negative Bacilli. (1988) N.Z.J. Med. Lab. Technol., 8-12.
4. CULLEN K.C., KLOOSTERMAN R.E., SHALIS P.J., PIERSON C.L.  
Comparison of Three Commercial Systems for the Identification of Glucose Non-fermenting Gram-Negative Rods. (1989) ASM Annual Meeting - Poster N° C26.
5. DANCE D.A.B., WUTHIEKANUN V., NAIGOWIT P., WHITE N.J.  
Identification of *Pseudomonas pseudomallei* in Clinical Practice: Use of Simple Screening Tests and API 20 NE. (1989) J. Clin. Pathol. 42, 645-648.
6. FRENEY J., GARONNAT D., BOUVARD V., FLEURETTE J.  
Comparaison de Deux Systèmes d'Identification de Bacilles Gram Négatifs Non Fermentants et de Bacilles Fermentants Oxydase Positive. (1984) Ann. Biol. Clin. 42, 337-341.
7. GEISS H.K., PIOTROWSKI H.D., HINGST V.  
Evaluation of API 20 NE in Routine Diagnostics of Nonfermenting Gram-Negative Rod-Shaped Bacteria. (1985) Zbl. Bakt. Mikrobiol. Hyg. A-Med. 259, 1, 35-42.
8. KISKA D.L., KERR A., JONES M.C., CARACCILO J.A., ESKRIDGE B., JORDAN M., MILLER S., HUGHES D., KING N., GILLIGAN P.H.  
Accuracy of Four Commercial Systems for Identification of *Burkholderia cepacia* and Other Gram-Negative Nonfermenting Bacilli recovered from Patients with Cystic Fibrosis. (1996) J. Clin. Microbiol. 34, 4, 886-891.
9. LAMPE A.S., VAN DER REIJDEN T.J.K.  
Evaluation of Commercial Test Systems for the Identification of Nonfermenters. (1984) Eur. J. Clin. Microbiol. 3, 301-305.
10. MARTIN R., SIAVOSHI F., McDOUGAL D.L.  
Comparison of Rapid NFT System and Conventional Methods for Identification of Nonsaccharolytic Gram-Negative Bacteria. (1986) J. Clin. Microbiol. 24, 1089-1092.
11. OVERMAN T.L., KESSLER J.F., SEABOLT J.P.  
Comparison of API 20 E, API Rapid E and API Rapid NFT for Identification of Members of the Family *Vibrionaceae*. (1985) J. Clin. Microbiol. 22, 778-781.
12. PALMIERI M.J., CARITO S.L., MEYER R.F.  
Comparison of Rapid NFT and API 20 E with Conventional Methods for Identification of Gram-Negative Nonfermentative Bacilli from Pharmaceuticals and Cosmetics. (1988) Applied and Environmental Microbiol. 54, 2838-2841.
13. PELADAN F., MONTEIL H.  
Identification of *Pseudomonas*, *Flavobacterium* and *Alcaligenes* with the API 20 NE System. (1988) Path. Biol., 36, 187-192.
14. SOGAARD P., GAHRN-HANSEN B., HUI-PING Z., FREDERIKSEN W.  
An Investigation of Three Commercial Methods for Rapid Identification of Non-Enteric Gram-Negative Rods. (1986) Acta Path. Microbiol. Immunol. Scand. Sect. B, 94, 357-363.
15. TOWNER K.J., CHOPADE B.A.  
Biotyping of *Acinetobacter calcoaceticus* using the API 20 NE System. (1987) Journal of Hospital Infection. 10, 145-151.
16. VON GRAEVENITZ A., ZOLLINGER-ITEN J.  
Evaluation of Pertinent Parameters of a New Identification System for Non-Enteric Gram-Negative Rods. (1985) Eur. J. Clin. Microbiol. 4, 108-112.
17. Clinical and Laboratory Standards Institute, M50-A Quality Control for Commercial Microbial Identification Systems; Approved Guideline, Vol. 28 n° 23.

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BIOMERIEUX, the BIOMERIEUX logo, ATB, API, APIWEB and ATB NEW are used, pending, and/or registered trademarks belonging to bioMérieux, or one of its subsidiaries, or one of its companies.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

Any other name or trademark is the property of its respective owner.

For users in the European Union (Regulation (EU) 2017/746) and in countries with similar requirements: Should a serious incident occur during the use of this device or as a result of its use, please report it to the manufacturer and/or their authorized representative as well as to your national authority

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